



## PHOTOFRIN® (porfimer sodium) for Injection

PHOTOFRIN® (porfimer sodium) for Injection is available exclusively through CuraScript, Inc.

## OPTIGUIDE™ Fiber Optic

The OPTIGUIDE® optical diffuser fiber is a light delivery system consisting of a 400µm core silica optical fiber, a proximal SMA-type connector, and a distal light diffusing tip. The fiber is designed to distribute light energy uniformly over the specified length of the diffuser tip.

The OPTIGUIDE® series fiber is intended for single use only. All OPTIGUIDE® series fibers are intended for use with 630 PDT lasers.

CuraScript Item Number	Diffuser Length	Type	Tip Outer Diameter	List Price (Each)
245283	1.0cm	Rigid	1.7mm	\$895.00
245291	1.5cm	Rigid	1.7mm	\$895.00
245305	2.0cm	Rigid	1.7mm	\$895.00
245313	2.5cm	Rigid	1.7mm	\$895.00
245321	5.0cm	Rigid	1.7mm	\$895.00
255432	1.0cm	Flexible	1.2mm	\$895.00
427381	2.5cm	Flexible	1.2mm	\$895.00
255459	5.0cm	Flexible	1.2mm	\$895.00

PHOTOFRIN® and OPTIGUIDE® are registered trademarks of Concordia Laboratories Inc. Pinnacle Biologics™ and the logo of Pinnacle Biologics™ are trademarks of Pinnacle Biologics, Inc. PHOTOFRIN® is distributed in the United States by Pinnacle Biologics, Inc.

US/PHO/NP/0018 (02/2024)

Pinnacle Biologics, Inc. Bannockburn, IL 60015 All rights reserved.

## ORDERING INFORMATION

To order OPTIGUIDE® fibers or PHOTOFRIN®, please contact CuraScript, Inc.

**Phone:** (866) 268-5554

**Customer Service:** (877) 599-7748

**Fax:** (800) 862-6208

**Website:** [www.curascriptonline.com](http://www.curascriptonline.com)

## CuraScript, Inc. can be contacted:

### Monday - Friday

8:30 a.m. – 7:00 p.m. Eastern Time

(Orders will be delivered the next business day)

## Product Inquiries

Reports about safety, adverse events, side effects and product complaints can be made by calling

1-866-248-2039 or by emailing

[medicalinformation@advanzpharma.com](mailto:medicalinformation@advanzpharma.com)

